

## Product datasheet for **HP219166**

### **C9orf164 Human qPCR Primer Pair (NM\_182635)**

#### **Product data:**

|                          |   |
|--------------------------|---|
| <b>Product Type:</b>     | qPCR Primer Pairs   |
| <b>Gene ID:</b>          | 349236  |
| <b>Forward Sequence:</b> | TCTCCTCAGCAGGGAACAAGAC  |
| <b>Reverse Sequence:</b> | TTCCTCCAGCTCTGCATCATCC  |
| <b>ACCN:</b>             | <a href="#">NM_182635</a> , <a href="#">NM_182635.1</a>   |
| <b>Synonyms:</b>         | FLJ39737  |
| <b>Components:</b>       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.   |
| <b>Quality Control:</b>  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| <b>Storage:</b>          | Store at -20 °C.  |
| <b>Stability:</b>        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)