

Product datasheet for **HP219173**

ANKDD1A Human qPCR Primer Pair (NM_182703)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	348094
Forward Sequence:	GCCACTTACGAATCCTCCAGATC
Reverse Sequence:	GGTCCTCCATTATGAACGCCAG
ACCN:	<u>NM_182703</u> , <u>NM_182703.1</u> , <u>NM_182703.2</u> , <u>NM_182703.3</u> , <u>NM_182703.4</u> , <u>NM_182703.5</u> , <u>BC027350</u> , <u>BC042002</u> , <u>BC047571</u> , <u>BC047607</u> , <u>BC101273</u> , <u>BC101274</u> , <u>BC101275</u> , <u>BC101276</u> , <u>BC110610</u> , <u>BX116779</u>
UniProt ID:	<u>Q495B1</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)