

## Product datasheet for **HP219181**

### IBRDC2 (RNF144B) Human qPCR Primer Pair (NM\_182757)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	255488
Forward Sequence:	GCAAACCTCTGCCTGTGTGAGCA
Reverse Sequence:	CATCCTTCTCGGATTGCCAGCT
ACCN:	<u>NM_182757, NM_182757.1, NM_182757.2, NM_182757.3, BC063311</u>
UniProt ID:	<u>Q7Z419</u>
Synonyms:	bA528A10.3; IBRDC2; p53RFP; PIR2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)