

Product datasheet for **HP219196**

P4HA3 Human qPCR Primer Pair (NM_182904)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283208
Forward Sequence:	AGAGTCACTGGCTCTGCCATCA
Reverse Sequence:	ATCCTCGGAAGAGACTGACAGC
ACCN:	<u>NM_182904</u> , <u>NM_182904.1</u> , <u>NM_182904.3</u> , <u>NM_182904.4</u> , <u>BC040582</u> , <u>BC089446</u> , <u>BC117333</u> , <u>BC126170</u> , <u>NM_182904.5</u>
UniProt ID:	<u>Q7Z4N8</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)