

Product datasheet for **HP219199**

RIMS4 Human qPCR Primer Pair (NM_182970)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	140730
Forward Sequence:	ACGGTCAGTTGGAGGTGGACAT
Reverse Sequence:	GATGCCATTCTCTAGCAGGTAGG
ACCN:	<u>NM_182970</u> , <u>NM_182970.2</u> , <u>NM_182970.3</u> , <u>BC036612</u> , <u>BC104190</u> , <u>BC104191</u> , <u>BC112985</u> , <u>BC160058</u> , <u>BU167487</u> , <u>BX648947</u> , <u>NM_182970.4</u>
UniProt ID:	<u>Q9H426</u>
Synonyms:	C20orf190; RIM-4; RIM 4; RIM4; RIM4-gamma; RIM4gamma
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)