

Product datasheet for **HP219216**

NBPF11 Human qPCR Primer Pair (NM_183372)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	200030
Forward Sequence:	GCATGGCTGTTGACATAGGCAG
Reverse Sequence:	CTATCCAGTGAGTCCTGCAAGAC
ACCN:	NM_183372 , NM_183372.1 , NM_183372.2 , NM_183372.3 , NM_183372.4 , NM_183372.5 , BC136292 , BI552657 , BM802549 , BP873295 , BX648497
UniProt ID:	Q86T75
Synonyms:	NBPF24
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)