

## Product datasheet for **HP219219**

### ESSPL (PRSS48) Human qPCR Primer Pair (NM\_183375)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	345062
Forward Sequence:	GTCACCTTCACTTCTGCCATCC
Reverse Sequence:	GCTTCTGAAGGGCAGAATGGT
ACCN:	<u>NM_183375, NM_183375.1, NM_183375.2, BC121109, BC121110, BN000134, NM_183375.4</u>
UniProt ID:	<u>Q7RTY5</u>
Synonyms:	ESSPL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)