

Product datasheet for **HP219225**

H1-9P Human qPCR Primer Pair (NM_194072)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	373861
Forward Sequence:	TTCCACCACGGGCTACGACAT
Reverse Sequence:	TCTGCTGTCGTTGCCTCTTGAC
ACCN:	NM_194072 , NM_194072.1 , BC033456 , BC048250 , BC109047 , BC109048
Synonyms:	histone H1-like protein in spermatids 1; OTTMUSP00000001456; TISP64
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)