

Product datasheet for **HP219231**

MORN2 Human qPCR Primer Pair (NM_194270)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	378464
Forward Sequence:	GGACTTACACATTCCCAAATGGG
Reverse Sequence:	TACCACTCCATTCTAGTCCTTGG
ACCN:	<u>NM_194270</u> , <u>NM_194270.1</u>
Synonyms:	MGC107136; Mopt; MOPT, BLOCK27, MGC126130, MGC126131, MGC126132; MORN repeat containing 2; protein containing single MORN motif in testis
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)