

## Product datasheet for **HP219254**

### KIF24 Human qPCR Primer Pair (NM\_194313)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	347240
Forward Sequence:	ATCCGAGCACTGGATCAGGAAC
Reverse Sequence:	GCTTGGTGAGATGTTGGCGATC
ACCN:	<a href="#">NM_194313</a> , <a href="#">NM_194313.1</a> , <a href="#">NM_194313.2</a> , <a href="#">NM_194313.3</a> , <a href="#">BC048311</a> , <a href="#">BC110502</a> , <a href="#">BC110503</a> , <a href="#">BE856625</a> , <a href="#">BF895212</a> , <a href="#">BP871643</a> , <a href="#">BU660279</a> , <a href="#">NM_194313.4</a>
UniProt ID:	<a href="#">Q5T7B8</a>
Synonyms:	bA571F15.4; C9orf48
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)