

Product datasheet for **HP219289**

GPR103 (QRFPR) Human qPCR Primer Pair (NM_198179)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	84109
Forward Sequence:	CCTGCTCATCACCTTCTTCTGC
Reverse Sequence:	GTCACAACAGCGGTAGACTGGA
ACCN:	<u>NM_198179</u> , <u>NM_198179.1</u> , <u>BC128132</u> , <u>BC128133</u> , <u>BG169612</u> , <u>NM_198179.3</u>
UniProt ID:	<u>Q96P65</u>
Synonyms:	AQ27; GPR103; SP9155
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)