

Product datasheet for **HP219342**

VSTM1 Human qPCR Primer Pair (NM_198481)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284415
Forward Sequence:	TCATCTGAGGAATCCACCAAGAG
Reverse Sequence:	GCTCAGCATAGGTCCTCCTTG
ACCN:	<u>NM_198481</u> , <u>NM_198481.1</u> , <u>NM_198481.2</u> , <u>NM_198481.3</u> , <u>BC100942</u> , <u>BC100943</u> , <u>BG189223</u> , <u>BQ446114</u>
UniProt ID:	<u>Q6UX27</u>
Synonyms:	SIRL-1; SIRL1; UNQ3033
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)