

## Product datasheet for **HP219347**

### FAM83H Human qPCR Primer Pair (NM\_198488)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	286077
Forward Sequence:	GAAGTTCCTGCTGGTGGACTGT
Reverse Sequence:	AGACCAGCTCTCCTTGGAACAC
ACCN:	<a href="#">NM_198488</a> , <a href="#">NM_198488.1</a> , <a href="#">NM_198488.2</a> , <a href="#">NM_198488.3</a> , <a href="#">BC007264</a> , <a href="#">BC033256</a> , <a href="#">BE856638</a> , <a href="#">BM982083</a> , <a href="#">NM_198488.4</a>
UniProt ID:	<a href="#">Q6ZRV2</a>
Synonyms:	AI3; AI3A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)