

Product datasheet for **HP219355**

VWA2 Human qPCR Primer Pair (NM_198496)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	340706
Forward Sequence:	ATGTGGCACTGCCATCCAAGCA
Reverse Sequence:	ATCCTCCACCTGCTCAGCCAAC
ACCN:	<u>NM_198496</u> , <u>NM_198496.1</u> , <u>BC128588</u> , <u>BC143380</u> , <u>BF355787</u> , <u>BF754676</u> , <u>BF802655</u> , <u>BF849816</u>
UniProt ID:	<u>Q5GFL6</u>
Synonyms:	AMACO; CCSP-2; NET42
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)