

Product datasheet for **HP219370**

GALNTL4 (GALNT18) Human qPCR Primer Pair (NM_198516)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	374378
Forward Sequence:	GGTGGATGACAACAGCAGTAACG
Reverse Sequence:	GCTTGCTGTGACGCACGACTTT
ACCN:	<u>NM_198516</u> , <u>NM_198516.1</u> , <u>NM_198516.2</u> , <u>BC060864</u> , <u>BC060864.1</u> , <u>BC037341</u> , <u>NM_198516.3</u>
UniProt ID:	<u>Q6P9A2</u>
Synonyms:	GalNAc-T15; GalNAc-T18; GALNACT18; GALNT15; GALNTL4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)