

## Product datasheet for **HP219386**

### ZNF568 Human qPCR Primer Pair (NM\_198539)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	374900
Forward Sequence:	AAGAGTCTGCCCTTTCCGAGGA
Reverse Sequence:	GCAGGTTTCATTTGCTCCCACTC
ACCN:	<u><a href="#">NM_198539</a></u> , <u><a href="#">NM_198539.1</a></u> , <u><a href="#">NM_198539.2</a></u> , <u><a href="#">NM_198539.3</a></u> , <u><a href="#">BC016334</a></u> , <u><a href="#">BC031218</a></u> , <u><a href="#">BC041927</a></u> , <u><a href="#">BC108696</a></u> , <u><a href="#">BC111777</a></u> , <u><a href="#">BC131547</a></u> , <u><a href="#">BX640681</a></u>
UniProt ID:	<u><a href="#">Q3ZCX4</a></u>
Synonyms:	ZFP568
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)