

Product datasheet for **HP219439**

GPR133 (ADGRD1) Human qPCR Primer Pair (NM_198827)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283383
Forward Sequence:	CAACCGAGTCTTCGTGTACTGC
Reverse Sequence:	GAGGATGGCAAAGTTGGTGAGG
ACCN:	<u>NM_198827</u> , <u>NM_198827.1</u> , <u>NM_198827.2</u> , <u>NM_198827.3</u> , <u>NM_198827.4</u> , <u>BC047951</u> , <u>BC112307</u> , <u>BC112309</u> , <u>BC143775</u> , <u>BQ028954</u> , <u>BX649116</u> , <u>NM_198827.5</u>
UniProt ID:	<u>Q6QNK2</u>
Synonyms:	GPR133; PGR25
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)