

Product datasheet for **HP219453**

LCN6 Human qPCR Primer Pair (NM_198946)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	158062
Forward Sequence:	GGCTTTGCCATGGAGAAGGACA
Reverse Sequence:	TCAGGTCCATGACACTCTGGTC
ACCN:	NM_198946 , NM_198946.1 , NM_198946.2 , BC062746 , BC062746.1 , BC040937
UniProt ID:	P62502
Synonyms:	hLcn5; LCN5; UNQ643
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)