

Product datasheet for **HP219457**

KCTD1 Human qPCR Primer Pair (NM_198991)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284252
Forward Sequence:	ACGCACGTCATCAGGTTCCAC
Reverse Sequence:	CGTATTCGCTGAACTGGGACGA
ACCN:	<u>NM_198991</u> , <u>NM_198991.1</u> , <u>NM_198991.2</u> , <u>NM_198991.3</u> , <u>BC063652</u> , <u>BC042371</u> , <u>BI667790</u> , <u>BI915955</u> , <u>BQ129398</u> , <u>BU153387</u> , <u>BX366041</u>
UniProt ID:	<u>Q719H9</u>
Synonyms:	C18orf5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)