

Product datasheet for **HP219461**

LIPI Human qPCR Primer Pair (NM_198996)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	149998
Forward Sequence:	TACAAGACTAGCTTATGTGTGGAC
Reverse Sequence:	CAGTGGTCCTAAGAGGTCTTCC
ACCN:	NM_198996 , NM_198996.2 , NM_198996.3 , BC140336 , BC148691 , BI460152
UniProt ID:	Q6XZB0
Synonyms:	CT17; LPDL; mPA-PLA1 beta; PLA1C; PRED5
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)