

Product datasheet for **HP219501**

LRR37A3 Human qPCR Primer Pair (NM_199340)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	374819
Forward Sequence:	AGAGCAGAGGCTCCGCCAAGT
Reverse Sequence:	CTCGGTCCAACGTATGCTTTCC
ACCN:	<u>NM_199340</u> , <u>NM_199340.1</u> , <u>NM_199340.2</u> , <u>NM_199340.3</u> , <u>BC043145</u> , <u>BC043145.1</u> , <u>BC048966</u> , <u>BC113643</u> , <u>BC117274</u> , <u>BI560889</u> , <u>NM_199340.4</u>
UniProt ID:	<u>O60309</u>
Synonyms:	LRR37; LRR37A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)