

## Product datasheet for **HP219520**

### SLC25A35 Human qPCR Primer Pair (NM\_201520)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	399512
Forward Sequence:	TGCACAAGGCTCTACAACCAGC
Reverse Sequence:	CCAGTAAAGGCGCTGAAGCACT
ACCN:	<u><a href="#">NM_201520</a></u> , <u><a href="#">NM_201520.1</a></u> , <u><a href="#">NM_201520.2</a></u> , <u><a href="#">BC101329</a></u> , <u><a href="#">BC041597</a></u> , <u><a href="#">BC052233</a></u> , <u><a href="#">BC101330</a></u> , <u><a href="#">BC105995</a></u> , <u><a href="#">BM987725</a></u> , <u><a href="#">NM_201520.3</a></u>
UniProt ID:	<u><a href="#">Q3KQZ1</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)