

Product datasheet for **HP219567**

IER5L Human qPCR Primer Pair (NM_203434)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	389792
Forward Sequence:	AGAACGGCTACTTGCACCAGGA
Reverse Sequence:	GCCAGGGTAATACTTGCCTTG
ACCN:	<u>NM_203434</u> , <u>NM_203434.1</u> , <u>NM_203434.2</u> , <u>BC013070</u> , <u>BC064028</u> , <u>BM763373</u> , <u>BQ574938</u> , <u>BU676697</u> , <u>BX458317</u>
UniProt ID:	<u>Q5T953</u>
Synonyms:	bA247A12.2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)