

Product datasheet for **HP219570**

GLIDR Human qPCR Primer Pair (NM_203448)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	389741
Forward Sequence:	ACTTCAGCACCACTGGCGGTG
Reverse Sequence:	CAGCGGAAAACGAGCAATGGGT
ACCN:	<u>NM_203448</u> , <u>NM_203448.1</u> , <u>BC000869</u> , <u>BC015151</u> , <u>BC019880</u> , <u>BC027338</u> , <u>BC040170</u>
Synonyms:	LINC01172; TCONS_00015562
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)