

Product datasheet for **HP219576**

NT2NL (NOTCH2NL) Human qPCR Primer Pair (NM_203458)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388677
Forward Sequence:	TGTCTCGACCTTGCTGAATGG
Reverse Sequence:	CATCGGTCCATTGGCACTCCTT
ACCN:	<u>NM_203458</u> , <u>NM_203458.1</u> , <u>NM_203458.2</u> , <u>NM_203458.3</u> , <u>NM_203458.4</u> , <u>BC065219</u> , <u>BC065219.1</u> , <u>BC010154</u> , <u>BC019835</u> , <u>BC084572</u> , <u>BC090943</u> , <u>BM792637</u> , <u>BX537434</u>
UniProt ID:	<u>Q7Z3S9</u>
Synonyms:	N2N; NOTCH2NL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)