

Product datasheet for **HP219585**

QIL1 (C19orf70) Human qPCR Primer Pair (NM_205767)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	125988
Forward Sequence:	TACCAGTTCAGCCAGTACGTGTG
Reverse Sequence:	CTGCATTCCAGGAGTCACGGAT
ACCN:	<u>NM_205767</u> , <u>NM_205767.1</u> , <u>NM_205767.2</u> , <u>BC009557</u> , <u>BC009557.1</u> , <u>BC042386</u> , <u>BC052640</u> , <u>BI827162</u> , <u>BM563199</u> , <u>NM_205767.3</u>
UniProt ID:	<u>Q5XKP0</u>
Synonyms:	C19orf70; MIC13; P117; QIL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)