

Product datasheet for **HP219594**

SPACA5 Human qPCR Primer Pair (NM_205856)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	389852
Forward Sequence:	CACCAAGAACCTCTGCCACATG
Reverse Sequence:	GAAAGCCCATTCTGTGAGGACAC
ACCN:	<u>NM_205856, NM_205856.1, NM_205856.2, NM_205856.3, BX279755</u>
UniProt ID:	<u>Q96QH8</u>
Synonyms:	LYC5; LYZL5; PNP6288; SLLP-X; SLLP2; UNQ6288
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)