

Product datasheet for **HP219601**

CCDC18 Human qPCR Primer Pair (NM_206886)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	343099
Forward Sequence:	CAAGGAAGACCGTTGCATTGGC
Reverse Sequence:	GCAGTCAGGTTTGCATGAAGCTT
ACCN:	<u>NM_206886</u> , <u>NM_206886.1</u> , <u>NM_206886.2</u> , <u>NM_206886.3</u> , <u>NM_206886.4</u> , <u>BC056667</u> , <u>BC152762</u> , <u>BC156745</u> , <u>BM149875</u> , <u>BX099893</u> , <u>BX538286</u>
UniProt ID:	<u>Q5T9S5</u>
Synonyms:	NY-SAR-41
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)