

Product datasheet for **HP219621**

C17orf87 (SCIMP) Human qPCR Primer Pair (NM_207103)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388325
Forward Sequence:	GTAAGTGGCAGCTTAGACGAGG
Reverse Sequence:	AGAGGCGGTAATTGAACTGGCG
ACCN:	<u>BC029580</u> , <u>NM_207103</u> , <u>NM_207103.1</u> , <u>NM_207103.2</u> , <u>NM_207103.3</u> , <u>BC130554</u> , <u>BC127925</u> , <u>BC130556</u> , <u>BC144395</u>
UniProt ID:	<u>Q6UWF3</u>
Synonyms:	C17orf87; UNQ5783
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)