

## Product datasheet for **HP219669**

### FRG1BP Human qPCR Primer Pair (NM\_207350)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	284802
Forward Sequence:	GACCTGGAACAGGATAGACG
Reverse Sequence:	AAGCTGGTTTCTTCATCTTCTCT
ACCN:	<a href="#">NM_207350</a> , <a href="#">NM_207350.1</a> , <a href="#">BC062724</a> , <a href="#">BC062743</a> , <a href="#">BC095491</a> , <a href="#">BU675133</a>
Synonyms:	bA348114.2; C20orf80; FRG1B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)