

Product datasheet for **HP219675**

LINC01140 Human qPCR Primer Pair (NM_207357)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	339524
Forward Sequence:	CCAGGCTTCTTCCACTTCCGAA
Reverse Sequence:	CTCCTTGTCTCCGCTCACTTTC
ACCN:	<u>NM_207357</u> , <u>BC030523</u> , <u>BC041867</u> , <u>BC126240</u> , <u>BC132861</u> , <u>BI818944</u> , <u>BM679066</u> , <u>BM722117</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)