

## Product datasheet for **HP219692**

### OCC1 (C12orf75) Human qPCR Primer Pair (NM\_207376)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	387882
Forward Sequence:	CCAAGGAAATTCTAAACAGTCACC
Reverse Sequence:	CCACAGAAAGTTCATAAACAAGTG
ACCN:	<a href="#">NM_207376</a> , <a href="#">NM_207376.1</a> , <a href="#">BC013920</a> , <a href="#">BC045812</a> , <a href="#">BC061920</a> , <a href="#">BC067897</a> , <a href="#">BC126980</a> , <a href="#">BC127012</a> , <a href="#">BC131552</a> , <a href="#">BQ021766</a>
UniProt ID:	<a href="#">Q8TAD7</a>
Synonyms:	AGD3; OCC-1; OCC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)