

## Product datasheet for **HP219704**

### CLEC17A Human qPCR Primer Pair (NM\_207390)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388512
Forward Sequence:	CTTGAGCCTTCTCCATTGCAGC
Reverse Sequence:	GAAGTCACCACCAGCAGACACA
ACCN:	<u><a href="#">NM_207390</a>, <a href="#">NM_207390.1</a>, <a href="#">NM_207390.2</a>, <a href="#">NM_207390.3</a>, <a href="#">BC140848</a>, <a href="#">BC144664</a>, <a href="#">BC144665</a></u>
UniProt ID:	<u><a href="#">Q6ZS10</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)