

Product datasheet for **HP219705**

RGS9BP Human qPCR Primer Pair (NM_207391)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	388531
Forward Sequence:	CTCCTCGGTCGTGCCTTGCA
Reverse Sequence:	CTTCGCCACGCACACGGCTAG
ACCN:	<u>NM_207391</u> , <u>NM_207391.1</u> , <u>NM_207391.2</u> , <u>BC172343</u> , <u>NM_207391.3</u>
UniProt ID:	<u>Q6ZS82</u>
Synonyms:	PERRS; R9AP; RGS9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)