

Product datasheet for **HP219731**

PADI6 Human qPCR Primer Pair (NM_207421)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	353238
Forward Sequence:	GTGGCATCTGTCTATGAGGACC
Reverse Sequence:	CTGAGGTGTGTCGAGGATCAAG
ACCN:	<u>NM_207421</u> , <u>NM_207421.1</u> , <u>NM_207421.2</u> , <u>NM_207421.3</u> , <u>NM_207421.4</u> , <u>BC148641</u> , <u>BC153197</u>
UniProt ID:	<u>Q6TGC4</u>
Synonyms:	hPADVI; PREMBL2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)