

Product datasheet for **HP219802**

HERV-FRD (ERVFRD-1) Human qPCR Primer Pair (NM_207582)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	405754
Forward Sequence:	ACCAGCTACCTGGGCATATCAG
Reverse Sequence:	GAGGCATTGGTGAATCGACTGG
ACCN:	NM_207582 , NM_207582.1 , NM_207582.2 , BC068585 , BC068585.1 , NM_207582.3
UniProt ID:	P60508
Synonyms:	envFRD; ERVFRDE1; GLLL6191; HERV-FRD; HERV-W/FRD; UNQ6191
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)