

Product datasheet for **HP219860**

RPL23AP7 Human qPCR Primer Pair (NR_000029)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	118433
Forward Sequence:	AGGCAGTGTGAAAGGTGTCCG
Reverse Sequence:	GGTGCTCTCTGAGGCTATCTG
ACCN:	<u>NR_000029</u> , <u>NR_000029.1</u> , <u>BC000596</u> , <u>BC105600</u> , <u>BE464804</u> , <u>BG335018</u> , <u>BI464848</u> , <u>BM661697</u> , <u>BX643196</u>
Synonyms:	bA395L14.9; RPL23AL1; RPL23A_6_267
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)