

## Product datasheet for **HP219874**

### CYP2B7P Human qPCR Primer Pair (NR\_001278)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1556
Forward Sequence:	TCATCGACACCTACCTGCTCCA
Reverse Sequence:	GTGGTCTCAGTGCCAGCAAAGA
ACCN:	<a href="#">NR_001278</a> , <a href="#">NR_001278.1</a> , <a href="#">BC041174</a>
Synonyms:	CYP2B; CYP2B7; CYP2B7P1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)