

## Product datasheet for **HP219946**

### ANXA2P1 Human qPCR Primer Pair (NR\_001562)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	303
Forward Sequence:	CGGACACATCTGGTGACTTCTG
Reverse Sequence:	TCCTCTTCACCCAGCGTCATA
ACCN:	<u>NR_001562</u>
Synonyms:	ANX2L1; ANX2P1; LPC2A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)