

Product datasheet for **HP220519**

PVT1 Human qPCR Primer Pair (NR_003367)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5820
Forward Sequence:	TTACAGGCGTGTGCCACAAAGC
Reverse Sequence:	GCCTGTAATCCCAGCACGTTGA
ACCN:	NR_003367 , BC033263 , BC041065 , BM821521 , BX280136 , BX648921
Synonyms:	LINC00079; MIR1204HG; NCRNA00079; onco-lncRNA-100
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)