

## Product datasheet for **HP220521**

### RRN3P1 Human qPCR Primer Pair (NR\_003370)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	730092
Forward Sequence:	GACACAAGCAGCTTTTGAGCGG
Reverse Sequence:	AACCACTGAGGGCAGGCAAATC
ACCN:	<u><a href="#">BC009198</a></u> , <u><a href="#">NR_003370</a></u> , <u><a href="#">BC006441</a></u> , <u><a href="#">BC068999</a></u> , <u><a href="#">BC070157</a></u> , <u><a href="#">BC112115</a></u> , <u><a href="#">BC112117</a></u> , <u><a href="#">BC143685</a></u> , <u><a href="#">BX101439</a></u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)