

Product datasheet for **HP220628**

HTR5BP Human qPCR Primer Pair (NR_003947)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	645694
Forward Sequence:	GATCTACGAGGCGCCAAGTTT
Reverse Sequence:	CCTCCTTCTCTAGCAACCGCTT
ACCN:	<u>NR_003947</u>
Synonyms:	5-Ht5b; 5-HT5B, GPR134; 5-hydroxytryptamine (serotonin) receptor 5B; MGC141198; OTTMUSP00000023087; protein-tyrosine-phosphatase (EC 3.1.3.48) Ly-5 precursor (B-cell variant)
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)