

Product datasheet for **HP220705**

UCA1 Human qPCR Primer Pair (NR_015379)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	652995
Forward Sequence:	CACTCTTTGCCAGCCTCAGCTT
Reverse Sequence:	AGGTGTGAGTGGCGGTCTGAAT
ACCN:	NR_015379 , BC005351 , BC017567 , BC022824 , BC055082 , BC070340 , BC107722
Synonyms:	CUDR; LINC00178; NCRNA00178; onco-lncRNA-36; UCAT1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)