

Product datasheet for **HP221775**

BOLA2 Human qPCR Primer Pair (AF271775)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	552900
Forward Sequence:	GAGGTTGTGCTGTAGCCCATG
Reverse Sequence:	CCTCAGGAAGTTCAGACTGTG
ACCN:	<u>BC022832</u> , <u>BC104642</u> , <u>BC144688</u> , <u>BQ639539</u> , <u>BU590412</u>
UniProt ID:	<u>Q9H3K6</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)