

Product datasheet for **HP222105**

MEG3 Human qPCR Primer Pair (AK055725)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	55384
Forward Sequence:	TTTGTGCCCAAGGCTCCTGGA
Reverse Sequence:	AGGGACTCAAGGAGCCAGTTA
ACCN:	BC023543 , BC036882 , BC051294 , BC062783 , BC092509 , BQ270557 , BQ613976
Synonyms:	FP504; GTL2; LINC00023; NCRNA00023; onco-lncRNA-83; prebp1; PRO0518; PRO2160
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)