

Product datasheet for **HP222130**

JKAMP Human qPCR Primer Pair (AK056145)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	51528
Forward Sequence:	TGATGCTGCTCCGACCTCTTCT
Reverse Sequence:	ACCAACTGCCTGAAGCACGGTT
ACCN:	<u>BC010359</u> , <u>BC010359.1</u> , <u>BC005198</u> , <u>BU928810</u> , <u>BX161501</u>
UniProt ID:	<u>Q9P055</u>
Synonyms:	C14orf100; CDA06; HSPC213; HSPC327; JAMP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)