

Product datasheet for **HP222773**

C4ORF28 (PACRGL) Human qPCR Primer Pair (AK098692)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	133015
Forward Sequence:	GGTGAACAGTCACGAGTGCCTT
Reverse Sequence:	GGAGGACATTCCTGTAATC
ACCN:	<u>BC023002</u> , <u>BC023002.2</u> , <u>BG741791</u> , <u>BM763154</u> , <u>BU674750</u> , <u>BX104331</u> , <u>BX344401</u>
UniProt ID:	<u>Q8N7B6</u>
Synonyms:	C4orf28
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)