

Product datasheet for **HP223219**

FAM89A Human qPCR Primer Pair (BC009447)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 375061 |
| Forward Sequence: | CTGCACTTACGCTCTGGAGAAC |
| Reverse Sequence: | GAGGAGACAGGCAGTGACAAGT |
| ACCN: | <u>BC009447</u> , <u>BC009447.1</u> |
| UniProt ID: | <u>Q96GI7</u> |
| Synonyms: | C1orf153; family with sequence similarity 89, member A; MGC15887; OTTHUMP00000036071; RP11-423F24.2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)