

## Product datasheet for **HP223223**

### **C9orf123 (TMEM261) Human qPCR Primer Pair (BC009510)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	90871
<b>Forward Sequence:</b>	CTGTTGAAGACCTGCTGGAGCT
<b>Reverse Sequence:</b>	GCCGATGACCATCTGCGTAATG
<b>ACCN:</b>	<a href="#">BC009510</a> , <a href="#">BC009510.2</a> , <a href="#">BI550091</a> , <a href="#">BQ446045</a>
<b>UniProt ID:</b>	<a href="#">Q96GE9</a>
<b>Synonyms:</b>	C9orf123; TMEM261
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)