

Product datasheet for **HP223257**

PRRG4 Human qPCR Primer Pair (BC010052)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	79056
Forward Sequence:	CATCCATGCTCTTCAGCCGTCT
Reverse Sequence:	TCCTGCATCCTCCACAGAAGGT
ACCN:	BC010052 , BC010052.2 , BC063393 , BC063393.1
UniProt ID:	Q9BZD6
Synonyms:	PRGP4; proline rich Gla (G-carboxyglutamic acid) 4 (transmembrane); TMG4; transmembrane gamma-carboxyglutamic acid protein 4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)